

# e-pill CADEX VibraPlus 8 Alarm Vibrating Watch

www.epill.com  
1-800-549-0095

## Quick Setup Guide

To switch between 24-Hour format and AM/PM format:  
Tap ST/STP when the current time is shown

### CLOCK – Set time & date (required)

1. The display should show the time & the day of the week. If not, tap **MODE** repeatedly until both are shown.
2. Press and hold **LIGHT** for 5 seconds until seconds digits in the box begin to flash. Then let go.
3. Tap **MODE** and the minutes will flash. Tap **ST/STP** to change the minute. Tap **MODE** to save.
4. The hours will now flash. Tap **ST/STP** to change the hour, note the AM/PM indicator. Tap **MODE** to save.
5. The date will now flash. Tap **ST/STP** to change the date. Tap **MODE** to save.
6. The month will now flash. Tap **ST/STP** to change the month. Tap **MODE** to save.
7. The year will now flash. Tap **ST/STP** to change the year. Tap **MODE** to save.
8. Tap **LIGHT** to confirm time & date.

**Note:** To temporarily display the date while in Time Mode, tap **RESET**.

### ALARMS - Set alarms in chronological order starting with the earliest (required)

1. Tap **MODE** repeatedly until Alarm Mode is reached. "A1" is shown in the box.
2. Set Alarm 1, press and hold **LIGHT** until hour flashes. Tap **ST/STP** to change the hour. Tap **MODE** to save.
3. The minutes will now flash. Tap **ST/STP** to change the minute. Tap **MODE** to save.
4. Tap **LIGHT** to confirm Alarm 1 (A1).
5. Tap **RESET** to select a different alarm (A1-A8) and repeat steps 2 through 4 to set it.
6. Set alarm notification type. Tap **ST/STP** once for sound, twice for vibrate or three-times for both.
7. Tap **MODE** repeatedly to go back to Time & Date.

**Note:** To silence the alarm at alarm time, tap any button. We recommend tapping **LIGHT**.

**Note:** All alarms repeat automatically at the set time(s) the next day.

### LOCK – Disable the buttons to prevent accidental reprogramming (optional)

1. The display should show the time & the day of the week. If not, tap **MODE** repeatedly until both are shown.
2. Press and hold **MODE** and **RESET** together at the same time.
3. "LOCK" will appear in the top left-hand corner of the display.
4. **LIGHT** button will still function but **ST/STP**, **MODE** and **RESET** buttons have now been locked and disabled.

**Unlock:** To unlock and enable all buttons hold **MODE** and **RESET** together at the same time. "LOCK" will now disappear.

### Answers to Frequently Asked Questions

1. Change the battery if the buttons stop working, the alarm is faint, or if the watch alarms/vibrates erratically.
2. Illumination: Tapping the **LIGHT** button only illuminates the digits (not the entire display). The function will only be noticeable when in a very dark room. At all other times, even at low light, the display will be readable.

# Battery Change Instructions for Jewelers

Please go to your local jeweler with these instructions

**Battery type CR2032**

1) Unscrew all 4 screws from backplate.



2) Remove the backplate. The sound maker is adhered to its inner face.


**Do Not Remove.**




3) Using a toothpick or small screwdriver, lift the battery up from the indicated opening.



4) Remove the sticker (with metal conductor) from the battery. Be careful not to fold the sticker.



5) Place the sticker and conductor on a fresh battery. Make sure the sticker is centered.



6) Place the new battery with the attached sticker back into the battery slot.



7) Make sure the battery is flush with the white frame, and the metal spring is uncovered.



8) Replace the backplate, be sure "CADEX" has the same orientation as the front face of the watch.




# e-pill® CADEX® VibraPlus™

8 Alarm Vibrating Watch  
Item: 932349 / 932444

Thank you for purchasing the e-pill CADEX VibraPlus.  
Your watch has the following functions:

1. Up to 8 daily auto-repeating alarms
2. Smart countdown timer
3. Notification types: Sound, vibrate, or both
4. Alarm duration: 15 seconds
5. Button lock
6. Light (illuminates digits only)



## Warranty


The inner components of this watch are warranted to be free from defects in materials and manufacturer's workmanship for one (1) year from the original date of purchase. This warranty does not cover the band, buckle, battery, case, buttons, or the quartz crystal. Do not wear this watch in the shower or swimming pool, it is **NOT WATERPROOF**. This warranty expressly excludes incidental and consequential damage caused by use, misuse, neglect, accident, unauthorized modification, or any defects resulting from service by anyone other than the manufacturer. During the warranty period, and upon proof of purchase, this watch will be repaired or replaced (with the same or a similar model) at the manufacturer's option. The e-pill CADEX VibraPlus watch is only a supplement to any memory system that you currently use to remind yourself to take your medications or perform medical tasks. The user maintains full responsibility to administer medications or perform medical tasks at the appropriate times. e-pill, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product, mathematical inaccuracy, or loss of stored data.

e-pill, LLC  
49 Walnut Street, Building 4  
Wellesley, MA 02481, USA  
www.epill.com

1-888-726-8805 (International: +1 781 239 2941)  
Written by e-pill Medication Reminders in Boston, MA  
Printed in China

## 4 Display Modes

Cycle through the 4 different display modes by tapping the **MODE** button.



## AM/PM Format or 24-Hour Format

To switch between 24-Hour format and AM/PM format:  
Tap **ST/STP** when the time is displayed (Time Mode).

## Time and Date

1. The display should show the time & the day of week. If not, tap **MODE** until both are shown.
2. Press and hold **LIGHT** for 5 seconds until the seconds digits in the box begin to flash. Then let go.
3. Tap **MODE** and the minutes will flash.
4. Tap **ST/STP** to change the minute. Tap **MODE** to save.
5. The hours will now flash. Tap **ST/STP** to change the hour, note the AM/PM indicator. Tap **MODE** to save.
6. The date will now flash. Tap **ST/STP** to change the date. Tap **MODE** to save.
7. The month will now flash. Tap **ST/STP** to change the month. Tap **MODE** to save.
8. The year will now flash. Tap **ST/STP** to change the year. Tap **MODE** to save.
9. Tap **LIGHT** to confirm time & date.

## Add Alarms

Set alarms in order starting with the earliest.

1. Tap **MODE** repeatedly until Alarm Mode is reached. "A1" will be shown in the box.
  2. Set Alarm 1, press and hold **LIGHT** until the hour begins to flash. Then let go.
  3. Tap **ST/STP** to change the hour. Tap **MODE** to save.
  4. The minutes will now flash. Tap **ST/STP** to change the minute. Tap **MODE** to save.
  5. Tap **LIGHT** to confirm Alarm 1 (A1).
  6. Tap **RESET** to select a different alarm (A1-A8) and repeat steps 2 through 4 to set it.
  7. Set alarm notification type. Tap **ST/STP** once for sound, twice for vibrate or 3 times for both. Sound: Vibrate:
  8. Tap **MODE** repeatedly to go back to Time & Date.
- Silence the alarm:** To silence the alarm at alarm time, tap any button. We recommend tapping **LIGHT**. All alarms repeat automatically at the set time(s) the next day.



### Delete an Alarm

1. Tap **MODE** repeatedly until Alarm Mode is reached. "A1" will be shown in the box.
2. Tap **RESET** repeatedly until the alarm you would like to delete is shown in the box.
3. Press and hold **RESET** until the corresponding alarm bell disappears from the middle of the display. Repeat step 2 and 3 for each alarm you would like to delete.
4. When done, press **MODE** until you return to Time Mode. Display should show time & day of week. Your alarm(s) have been erased.

### Smart Countdown Timer

The countdown timer only counts down within the time window that you set. Within these hours the watch will alarm every time the countdown reaches zero and automatically restart after alarming. To pause the countdown until the next day press **ST/STP** while in TR mode. To start a paused countdown press **ST/STP** while in TR mode.

**Before you begin:** When starting the countdown interval for the first time, the current time will need to be within the set time window.

1. Set the notification type. Tap **MODE** repeatedly until Alarm Mode is reached. "A1" will be shown in the box.
  2. Tap **ST/STP** once for sound, twice for vibrate or 3 times for both. Sound:  Vibrate: 
  3. Tap **MODE** until "TR" is shown above the line (TR Mode).
  4. Press and hold **LIGHT** until the starting hour begins to flash, this is when the interval will start counting down every day.
  5. Tap **ST/STP** to change the starting hour. Tap **MODE** to save.
  6. The ending hour will now flash, this is when the interval will stop counting down every day.
  7. Tap **ST/STP** to change the ending hour. Tap **MODE** to save.
  8. The countdown duration hours will now flash, this is how many hours the countdown will last.
  9. Tap **ST/STP** to change the duration hours. Tap **MODE** to save.
  10. The duration minutes will now flash, this is how many extra minutes the countdown will last.
  11. Tap **ST/STP** to change the duration minutes. Tap **MODE** to save.
  12. Tap **LIGHT** to confirm countdown interval settings.
  13. Tap **ST/STP** to start the countdown (TR will begin to flash).
- When the countdown reaches zero:** To silence the alarm and automatically restart the countdown timer, tap any button. We recommend tapping **LIGHT**.

**Permanently turn off:** To permanently turn off the countdown timer set both the starting and ending hour to 12am. Steps 3-7 above. Then tap **MODE** to return to time mode.

### Lock

Disable the buttons to prevent accidental reprogramming. The **LIGHT** button will still function but **ST/STP**, **MODE** and **RESET** buttons will be locked and disabled.

1. The display should show the time & the day of week. If not, tap **MODE** until both are shown.
  2. Press and hold **MODE** and **RESET** together at the same time.
  3. "LOCK" will appear in the top left-hand corner of the display.
- Unlock:** To unlock and enable all buttons hold **MODE** and **RESET** together at the same time. "LOCK" will now disappear.

### Stopwatch

1. Tap **MODE** until "ST" is shown above the line (ST Mode).
2. Tap **ST/STP** to start the stopwatch.
3. Tap **ST/STP** to stop the stopwatch.
4. Tap **RESET** to return to 0.

### Battery

Change the battery if the buttons stop working, the alarm is faint or if the watch alarms/vibrates erratically. When replacing the battery we strongly recommend that you take your e-pill watch to a jeweler with the battery change instructions that are found on page 2.

Your e-pill CADEX VibraPlus uses a CR2032 button cell battery. The CR2032 battery size is the most common watch battery in the world and can be found at any jeweler, grocery store or pharmacy.

### Care

DO NOT IMMERSE IN WATER. Avoid exposing the watch to temperature extremes. Avoid dropping the watch. Avoid exposing the watch to gasoline, cleaning solvents, adhesive agents, paint, or aerosol sprays.

### Caution

Dead battery, replacement of battery, malfunction or repairs will cause all data stored in the watch to be deleted. In no event will e-pill, LLC be liable for loss of data or damages arising from such loss. This watch will withstand splashes and rain but is not waterproof. DO NOT USE IN SHOWER, BATHTUB OR SWIMMING POOL.