

Medication Reminders

Operating Instructions

e-pill® MedGlider HOME

SKU: 843632 UPC: 837066007137

Thank you for buying
an e-pill medication reminder!

We suggest you read the following
instructions carefully and keep
them for future reference.

This easy-to-use e-pill reminder can remind
you to take your medications up to 4 times
per day. It has a flashing visual indicator and
a beeping alarm. This e-pill includes a unique
“missed pill” feature. The alarms automatically
reset every day at midnight.

e-pill® Medication Reminders

e-pill, LLC

49 Walnut Street, Bldg. 4

Wellesley, MA 02481, USA

1-800-549-0095 (+1 781 239 2941)

fax + 1 (781) 235-3252


www.epill.com

Version 140828

Setting up Your e-pill MedGlider HOME

If your device is not activated, please activate it now by removing the clear plastic tab from the back of your e-pill Medication reminder.

Setting Time and Date



Lift the clear
button cover
on the face of
the timer.

Hold down **SET** button until the time
shown on the display screen begins to flash.

TIME: Press the **HOUR** button to increase
hour (note AM / PM), Press **MIN** button to
increase minutes. Once desired time is
reached, press **SET**.


MONTH: The MONTH number will start
to flash on the bottom left corner of the
display screen. Press the **HOUR** button
to reach correct month # (1-12). Press **SET**.

DAY: The DAY of the month will start to
flash. Press the **HOUR** button to reach
the correct day (1-31). Press **SET**.

YEAR: The YEAR will start to flash. Press
the **HOUR** button to reach the correct
year (20xx). Press **SET**.

DONE: Your settings have been saved.

Setting Alarms



Month Day Year First
Dosage
Time Missed Pills

1. Press **MODE** button.
2. The 1 will start to flash on the display in
the lower right hand corner.
3. Press and hold down **SET** until - - : - -
starts to flash.
4. To change the hour, press the **HOUR** button.
Continue to press **HOUR** until the desired alarm
hour is reached. *Please observe the
AM/PM Indicator.
5. To change the minutes, press the **MIN** button.
Continue to press **MIN** until the desired alarm
minutes are reached.
6. Press **SET** to store the alarm time.
7. Press **MODE** and 2 (the second alarm)
will start to flash.
8. Repeat steps **3 - 6** for alarms 2 thru 4.

Daily Use

Turning Off Alarms

When the timer beeps with the red light flashing that means its time for you to take your medications.

Check to see which number (1-4) is flashing on the right bottom of the display screen.

If the number "1" is flashing, it is time to take your first dose of pills. Press the **STOP** button to stop the beeping & flashing light, and take the pills from pillbox "1".

The same should be done for alarm "2", "3", and "4".

Alarms will **automatically reset** at midnight each day and will repeat at the same time each day.

MISSED PILL INDICATOR

If "**STOP**" button is not pressed when the alarm is sounding, the alarm number will than be shown with an X meaning "**MISSED PILL**".

Any time pills are missed, you should follow the directions from your phycian or pharmacist before taking any medication (pill) out of sequence.

Erasing Alarms

To erase an alarm, simply press **MODE** until you reach the alarm you want to change (1-4) .

Once you have reached the correct alarm hold down **SET** until - - : - - appears.

Press **SET** twice and only the current alarm number (1-4) will blink on the screen.

Press **MODE** to reach the remaining alarms and time.

Replacing the Batteries



1. Slide open the battery door by pushing down on the arrow firmly with your thumb.

2. Insert 2 AAA size batteries, making sure the negative (-) and positive (+) terminals are placed correctly. We recommend using two alkaline batteries.

3. Put the battery compartment cover back and lock it into place.

Disclaimer

This e-pill Medication Reminder is only a supplement to any memory system that you currently use to remind yourself to take medications. The user maintains full responsibility to administer medication at the appropriate times. e-pill, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product. e-pill, LLC's liability shall in no circumstances exceed the amounts paid by the user for the timer.