

BATT: Set the type of battery in the dispenser for optimal battery life.

1. When in menu B, tap button 1 until 'BATT' is displayed.
2. Tap button 2 or 3 to set battery type to either:
  - a. 'ALC' for alkaline/non-rechargeable (most common).
  - b. 'NIMH' for rechargeable (uncommon).

SECURE: Auto-locking dispensing window. This setting refers to the dispensing window only, the window where the pills dispense from. The entire lid should always be kept locked.

ON: The dispensing window in the front of the lid only unlocks when the alarm is sounding. The window will lock again 2 minutes after dispensing.

OFF: The dispensing window in the front of the lid never locks.

1. When in menu B, tap button 1 until 'SECURE' is displayed.
2. Tap button 2 to turn on or button 3 to turn off.

HOURS: Use either a 12-hour clock (AM/PM) or 24-hour clock.

The dispenser's default uses the 24-hour clock.

1. When in menu B, tap button 1 until 'HOURS' is displayed.
2. Tap button 2 or 3 to set to:
  - a. 12 hour – 'AM' or 'PM' will be displayed next to the time.
  - b. 24 hour – 'AM' or 'PM' will not be displayed next to the time.

DLS: Set the dispenser to automatically adjust for daylight saving.

1. When in menu B, tap button 1 until 'DLS' is displayed.
2. Tap button 2 or 3 to set to:
  - a. 'OFF' to disable.
  - b. 'EU' for European Daylight Saving.
  - c. 'US' for United States Daylight Saving.

**DISPENSER START: DO NOT SET IF YOU ARE USING THIS PILL DISPENSER DAILY.** To delay the time at which the dispenser will start to operate. This function is useful if the dispenser is being filled, for example by a caregiver offsite and is not being used until a later date.

If this function is not needed, tap button 1 until you reach STOP.

If this function is needed, use the steps below to set the date and time of when the filled dispenser is getting delivered to the patient. When in menu B, tap button 1 until 'START' is displayed.

1. Set the Year (YY) Month (MM) Date (DD) and Hour (HH) that you wish the dispenser to start operating. The time (HH) should be in advance of the first dose to be dispensed for the user.
2. Tap button 2 (decrease) or 3 (increase) the settings.

DISPENSER STOP: It is possible to stop the dispenser's operation at a future date and time. This function is useful if two or more dispensers are being used by a user in sequence and ensures that only one dispenser is active at one time.

If this function is not needed, tap button 1 until you reach RECORD.

If this function is needed, use the steps below.

When in menu B, tap button 1 until 'STOP' is displayed.

1. Set the Year (YY), Month (MM), Date (DD) and Hour (HH) to when you want the dispenser to stop operating. The time (HH) should be set after the last dose to be dispensed for the user.
2. Tap button 2 (decrease) or 3 (increase) the settings. **If this function is not required set the year (YY) to a distance future year ex: YY38**

RECORD: Will show the time at which the previous 28 doses were dispensed

1. When in menu B, tap button 1 until 'RECORD' is displayed.
2. Tap button 2 (decrease) or 3 (increase) the previous dispensing times:
  - a. '1' is the latest dispense time which occurred.
  - b. '2' is the second latest dispense time which occurred.
  - c. '3' is the third latest dispense time which occurred etc.

SAVE and EXIT Menu B:

1. Press and hold button 1 for 5 seconds to save and exit Menu B.
2. The display will cycle between the current time and the current date.

## SETUP - STEP 2

### Filling the Tray

1. With the lid still open, remove the tray that is inside the dispenser. This makes filling the tray easier. Important: Never force the tray to rotate.
2. e-pill MedTime Pro comes with three different dosing templates. The template identifies the number of alarms per day and, very importantly, the time of day that you are filling the tray. The blank disc is used for 5 -