



SET UP INSTRUCTIONS EXTRA PENDANT Model #35915

Contents



Pendant



Rechargeable Li-Ion Special Pendant Battery - 1 piece



Belt Clip



Wrist / Walker Wheelchair strap



LogicMark v3.2

This Extra Pendant is to be used with the Freedom-Alert product model 35911. Up to 4 of these pendants (or a total of 4 pendants and Emergency Wall Communicators) can be taught to one base unit. For full instruction on the use and limitation of this product - refer to the FreedomAlert user manual.

Before programming this extra pendant to work with your FreedomAlert base unit, you must acknowledge that all of the Warnings, Limitations of Liability, Disclaimers, FCC information and Testing Instructions written in the FreedomAlert instruction manual pertain to the use of this product as well when used with the FreedomAlert product. By programming this extra pendant to your FreedomAlert product, you acknowledge that you have a copy of these instructions and have read them. If you do not have a copy, you can always get a copy off of our website or by calling our company and requesting one.

Installing the Battery

The Pendant requires one of the white Li-Ion (Lithium Ion) batteries at any given time. The pendant should be able to run up to 4 months on a fully charged battery. Simply unscrew the pendant battery cap and install the Pendant battery with the positive side down so that the negative side of the battery is touching the inside of the battery cap when the cap is on. You will hear a voice confirmation of the battery status.

Teaching the Extra Pendant to the FreedomAlert System.

On the Base Unit:

- Push and <u>release</u> the RED Code Learn button on the back of the base unit. You will hear a voice prompt "Pendant Learning"
- On the Pendant:
- 2. PUSH the Gray Battery Test button and the Blue Panic button at the same time-HOLD both buttons until you hear "Pendant Learning" from the pendant - then release both buttons on the pendant.
 - If "Base and Pendant out of Range" is heard from the pendant, STOP - wait 30 seconds and start at step #1 again.

- 3. When Pendant learning has been successful you will hear "All Systems are OK" and then it will announce the call setting that the base is set to Friends 911 Mode, Friends Only or 911 Only.
 - If "Base and Pendant out of Range" is heard from the Base or Pendant, STOP wait 30 seconds and start at step #1 again.

Notes:

- You have up to 45 seconds after pushing the RED Code Learn button on the base unit to press both buttons until you hear the second "Pendant Learning" announcement.
- Teaching additional pendants does not erase previously learned pendants.
- The base unit will automatically exit the pendant learning mode after several seconds or "Pendant Learning Failed".
 - If you hear the voice prompt "Base and Pendant Out of Range" an error has occurred, wait 30 seconds for the base to reset and repeat steps 1 & 2.
- A combined total of 4 pendants and Emergency Wall Communicators can be learned to a base unit. <u>If a pendant has been lost or</u> <u>replaced - ERASE all pendants and relearn</u>

the current pendants and any Emergency Wall Communicators.

Erasing all previously learned pendants:

Simply hold down the red **Code Learn** button at the back of the Base unit for 8 seconds. You will hear THREE voice prompts: "Erasing all pendants," "All pendants erased," and the system announcing the version and call setting of the unit. After that, you can program a new pendant to the base unit following the steps above as desired.

Battery Check / System Check:

The FreedomAlert has a built-in battery tester and method to check the whole system. You should make yourself familiar with these tests from the main FreedomAlert user's manual. If you do not have these important safety instructions, please get them from our website or by calling the company.

Water Resistance:

Please note that although the pendant is water resistant, it is not waterproof. Regular exposure to water, soaps, shampoos and conditioners can damage the water seals over time as well as affect the performance of the speaker – reducing the volume. To prolong the life of your pendant, we recommend laying the pendant just outside of the shower or bathtub – but well within reach should an accident happen in the shower or bathtub area.

If you have further questions on this product, call our TECH SUPPORT department at 703-934-7934 (option 2 for customer service).

Information The FCC Wants You To Know:

FCC Part 15 ID: TYD3X911

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Canada IC: 8471A-3X911 This Class B digital apparatus complies with Canadian ICES-003.

ONE YEAR LIMITED WARRANTY

LogicMark, LLC warrants to the original consumer/purchaser that this product shall be free of defects in material and workmanship under normal use and circumstances for a period of **one** (1) year from the date of original purchase for use. When the original consumer/purchaser returns the product pre-paid to LogicMark, LLC; 10106 Bluegrass Pkwy; Louisville KY, 40299, USA within the warranty period, and if the product is defective, LogicMark, LLC will at its option repair or replace such.

This warranty shall constitute the sole liability of LogicMark, LLC concerning the product. LogicMark, LLC expressly disclaims all other warranties INCLUDING, WITHOUT LIMITATION, THE WARRANTIES OF MERCHANT ABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NO PERSON, FIRM, OR CORPORA-TION IS AUTHORIZED TO ASSUME FOR LogicMark, LLC. ANY OTHER LIABILITY IN CONNECTION WITH THE SALE AND USE OF THE PRODUCT. LogicMark, LLC and agents and distributors will bear no liability whatsoever for incidental or consequential damages or charges of any kind.

Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above disclaimer regarding incidental or consequential damages may not apply to you.

This warranty is void if the product has been damaged or tampered with or if the product or any such parts have been opened. In all cases of damage during shipment, a claim must be filed with the shipping carrier and not with LogicMark, LLC.

This warranty gives you specific legal rights; you may also have other rights which vary from state to state.

OUT OF WARRANTY REPAIRS

If the warranty period has expired or if you are not the original owner of the product, LogicMark, LLC will at its option either (1) replace this product with a functionally similar (but not necessarily identical) refurbished product or (2) repair the original product and return it to the original consumer/purchaser after payment of repair/replacement charges has been received.

LogicMark

10106 Bluegrass Pkwy; Louisville KY, 40299 Tel: 1-703-934-7934 Toll Free: 1-800-519-2419 Fax: 1-703-934-7935 www.LogicMark.com

FreedomAlert is covered by the following patent: US # 8,369,821

©2013 LogicMark, LLC