
Setup Guide

e-pill® Monthly 5 Compartment Pillbox with 6 Alarm Pocket Reminder
Item: 960500




Up to Six Alarms


Easy to Set


Loud Alarms


Medical Snooze

TURN ON THE TIMER

1. Insert the 2 included AAA batteries.
2. The display will automatically turn on.

SET THE CLOCK

1. Open the clear lid.
2. Choose between 12-hour and 24-hour mode:
   - Press and hold the **STOP** button, then tap the **H** button.
     "**A**" will indicate **AM**, and "**P**" will indicate **PM**.
3. Set the hour:
   - Tap the **H** button until the desired hour is displayed.
4. Set the minutes:
   - Tap the **M** button until the desired minutes are displayed.
5. The time will be saved automatically and stop flashing 5 seconds after you release all buttons.

SET ALARMS (Up to 6 Alarms)

1. Set an alarm:

Push down the switch(es) for the desired alarm(s).

- Red switches set **AM** alarms (Midnight to 11 AM).
- Blue switches set **PM** alarms (Noon to 11 PM).

An alarm bell symbol and a number will appear on the screen for each set alarm.



2. Check Alarms:
   - Tap the **STOP** button to cycle through the set alarms.
   - "**A**" indicates **AM**, and "**P**" indicates **PM**.
3. Close the clear lid when finished.

ALARM REPEAT

If the alarm is not stopped within 5 minutes, it will repeat for 15 seconds one more time.

The corresponding alarm  bell will blink for 5 minutes or until the **STOP** button is pressed.

BATTERY INFORMATION

The timer uses 2 AAA batteries. Battery life is

approximately 1 year.

The battery compartment is located on the back of the timer.

STOP THE ALARM

1. When the alarm sounds, take your medications.

2. Tap the **STOP** button to silence the alarm.

Contact US

e-pill, LLC

49 Walnut Street., Bldg. 4

Wellesley, MA 02481, USA

Tel: 1-800-549-0095 | 1-781-239-2941

Fax: 1-781-235-3252

Email: sales@epill.com

Website: www.epill.com

Disclaimer

This e-pill medication reminder is only a supplement to any memory system that you currently use to remind yourself to take medications. The user maintains full responsibility to administer medication at the appropriate times. e-pill, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product.

e-pill, LLC's liability shall under no circumstance exceed the list price on epill.com.

Version 260211 | Printed in China