



ABPM 7100 HYPERTENSION MANAGEMENT SYSTEM (HMS) SPECIFICATIONS



SYSTEM REQUIREMENTS

COMPUTER:

- Processor: min. 1 GHz
- Memory: min. 2 GB RAM
- Hard disk space: min. 500 MB
- Resolution: min. 1024x768 pixels
- At least one free USB socket

OPERATING SYSTEM:

- Windows 10, Windows 11 (32-bit & 64-bit)
- macOS 13 Ventura (without Bluetooth support)

SOFTWARE:

- Installing the Software:
 - Java Runtime Environment

BlueSoleil drivers and software are not supported by IEM GmbH and may lead to complications with device communication.

Explore the entire Baxter cardiopulmonary line, including the Ambulatory Blood Pressure Monitor 7100 at hillrom.com/en/products/abpm-7100 or contact your Baxter Sales Representative.

Rx Only. For safe and proper use of the products mentioned herein, please refer to the appropriate Operator's Manual or Instructions for Use.

Baxter.com | Baxter International Inc.

Baxter is a trademark of Baxter International Inc.

US-FLC26-240001 (v1.0) 10/2024