



Nellcor™ OxiMax™ N-560™ Pulse Oximeter

Nellcor™
Puritan Bennett™
Airox™
Mallinckrodt™
DAR™
Shiley™
Sandman™

FEATURING LOSAT™ EXPANDED ACCURACY

The LoSat™ expanded accuracy feature of our OxiMax adhesive sensors assures you of a wide accuracy range (60% to 100% SpO₂) when used with the OxiMax N-560 pulse oximeter. This allows improved patient assessment at the more challenging lower saturation levels common in infants with congenital heart disease.

MONITOR WITH CONFIDENCE.

The Nellcor™ OxiMax™ N-560™ pulse oximeter from Covidien is a compact, affordable bedside monitor that delivers full-size performance. It incorporates the latest Nellcor signal processing technology for accurate, reliable SpO₂ and pulse rate measurements even in conditions of low perfusion and signal interference, including patient motion. This technology, along with the full line of versatile OxiMax sensors, allows clinicians to confidently and effectively monitor a wide range of patients.

Sleek, simple design. The OxiMax N-560 pulse oximeter is smaller than most other standalone pulse oximeters and has a simple, well-designed user interface that caregivers will be comfortable with. Save time and space with this easy-to-use OxiMax pulse oximeter.



OXIMAX N-560 PULSE OXIMETER FEATURES

- Monitors neonates through adults.
- Simple, intuitive operation and space-saving design.
- Nellcor's distinctive digital signal processing technology enhances performance during low perfusion, patient motion and other forms of signal interference.
- Variable pitch beep tone enables clinicians to hear point-by-point changes in SpO₃.
- Display changes from green to red during alarms.
- Eight-hour battery life.
- 24-hour trend memory.
- Optional paging system supplements bedside and central station alarms to reduce response time.
- Integrated Score analysis functionality stream lines in-hospital sleep screening.
- Patient trend data can be stored on a USB drive for archive and analysis.

TEAM UP WITH TELEMETRY

To enhance patient safety monitoring on your general care floor, you can team OxiMax N-560 pulse oximeters at the bedside with Covidien's cost-effective remote oximetry alarm systems or connect them to nurse call systems.

OXIMAX SENSOR ADVANTAGES

As part of Covidien's Nellcor™ OxiMax™ pulse oximetry system, the OxiMax™ N-560™ pulse oximeter works with the full line of OxiMax sensors, which help expand your patient monitoring options.

 OxiMax Max-Fast[™] forehead sensor offers improved performance when patients have low perfusion.

 OxiMax SoftCare[™] nonadhesive sensors provide a gentle, single-patient-use alternative to adhesive sensors for patients with fragile skin.

FEATURES AND SPECIFICATIONS

PERFORMANCE

MEASUREMENT RANGE

SpO₃: 1% to 100%

Pulse rate: 20 to 250 beats per minute (bpm)

Pulse amplitude: 0.03% to 20%

ACCURACY

Saturation (% $SpO_3 \pm 1 SD$)

Adult: 70% to 100% \pm 2 digits

60% to 80% \pm 3 digits

Neonate: 70% to 100% \pm 2 digits 60% to 80% \pm 4 digits

Low perfusion: 70% to 100% \pm 2 digits

Pulse rate

20 to 250 bpm \pm 3 digits

Low perfusion: 20 to 250 bpm \pm 3 digits

ELECTRICAL

INSTRUMENT

Power requirements: 100 to 240 VAC. 25 to 39 VA

Fuse rating: Qty 2, 2A, 250 volts,

slow-blow IEC (5x20 mm)

BATTERY

Type: NiMH

Battery capacity: Minimum of 8 hours using new,

fully charged battery with no

alarms, no serial data

ENVIRONMENTAL

OPERATING TEMPERATURE

10°C to 45°C Instrument:

Transport/Storage Temperature

(in shipping carton): -20°C to 70°C

Operating Humidity

15% to 95% noncondensing

Operating Altitude

-1,280 ft to 18,288 ft

PHYSICAL CHARACTERISTICS

WFIGHT

3.07 lb (1.39 kg)

SI7F

2.95 H x 9.06 W x 5.04 D (inches) (7.5 cm x 23.0 cm x 12.8 cm)

EOUIPMENT COMPLIANCE

• IEC 60601-1 • ISO 9919:2005

• CSA 601.1 • EN/IEC 60601-1-2

• UL 2601-1 (2nd ed.)

IPX1

EMISSIONS COMPLIANCE

• EN55011, Group 1, Class B

OUTPUT

- EIA-232 digital output from serial port
- Oxinet[™] Communications Output Protocol
- · Nurse call capability

DISPLAY/INDICATORS

High quality LED display changes from green to red

during alarm conditions.

Pulse amplitude indicator (10-segment LED).

Visual indicators: Pulse search, audible alarms silenced or off, interference indicator, battery charging, data in sensor and SatSeconds[™] alarm management clock

Audible and visual alarms for high/low saturation and pulse rate.

SatSeconds Alarm Management settings:

10, 25, 50 and 100, or OFF

Audible and visual warning indicators for low battery

and sensor off.

Audible sensor disconnect alarm.

OPTIONAL ACCESSORIES

- · GCX mounting adapter
- Pole mount clamp
- Interface cables
- Visual alarm indicator antenna

trademarks of a Covidien company. © 2009 COVIDIEN. All rights reserved. 09-PM-6420 MN 07608

Adult specifications are shown for the OxiMax™ MAX-A and MAX-N sensors with the OxiMax™ N-560™ pulse oximeter. Neonate specifications are shown for MAX-N sensors with the OxiMax[™] N-560[™] pulse oximeter. Saturation accuracy will vary by sensor type. Refer to Nellcor's[™] Sensor Accuracy Grid Card for details. The low perfusion specifications apply to monitor performance. Refer to the OxiMax[™] N-560[™] pulse oximeter operator's manual for complete description, instructions, warnings, cautions and specifications. Specifications are subject to change without notice

COVIDIEN, COVIDIEN with logo, Covidien logo and positive results for life are U.S. and/or internationally registered trademarks of Covidien AG. All other brands are



BOULDER, CO