

Chlorine, Free Portable Photometers

10

www.easeongear.com
[click to buy: www.hannameters.com](http://www.hannameters.com)


- CAL CHECK™
- User calibration
- Certified calibration and verification standards
- BEPS (Battery Error Prevention System)
- TIMER function
- Auto shut-off
- GLP Features

The HI 96701 meter measures the free chlorine (Cl_2) content in water samples in the 0.00 to 5.00 mg/L (ppm) range.

HI 96762 was specially developed to measure low concentrations of free chlorine in drinking water.

These meters use an exclusive positive-locking system to ensure that the cuvette is in the same place every time it is placed into the measurement cell.

Chlorine is widely used in making many everyday products, but most notably, it is used for producing safe drinking water the world over. Even the smallest water supplies are now usually chlorinated. It is also extensively used in the production of paper products, dyestuffs, textiles, petroleum products, medicines, antiseptics, insecticides, food, solvents, paints, plastics, and many other consumer products. Most of the chlorine produced is used in the manufacture of chlorinated compounds for sanitation, pulp bleaching, disinfectants, and textile processing.

Organic chemistry demands much from chlorine, both as an oxidizing agent and in substitution, since it often brings many desired properties in an organic compound when substituted for hydrogen, as in one form of synthetic rubber.

SPECIFICATIONS	HI 96701 Chlorine, Free	HI 96762 Chlorine, Free Low Range
Range	0.00 to 5.00 mg/L (ppm)	0.000 to 0.500 mg/L (ppm)
Resolution	0.01 mg/L from 0.00 to 3.50 mg/L (ppm); 0.10 mg/L above 3.50 mg/L (ppm)	0.001 mg/L (ppm)
Accuracy @ 25°C (77°F)	±0.03 mg/L ±3% of reading	±0.020 mg/L ±3% of reading
Light Source	tungsten lamp	
Light Detector	silicon photocell with narrow band interference filter @ 525 nm	
Power Supply	9V battery	
Auto-off	after ten minutes of non-use in measurement mode; after one hour of non-use in calibration mode; with last reading reminder	
Environment	0 to 50°C (32 to 122°F); RH max 95% non-condensing	
Dimensions	192 x 104 x 69 mm (7.6 x 4.1 x 2.7")	
Weight	360 g (12.7 oz.)	
Method	adaptation of the USEPA method 330.5 and Standard Method 4500-Cl G	

The reagents are in powder form and are supplied in packets. The amount of reagent is precisely dosed to ensure the maximum repeatability.

ORDERING INFORMATION

HI 96701 and HI 96762 are supplied with sample cuvettes (2) with caps, 9V battery and instruction manual.

CAL CHECK™ standards and testing reagents sold separately

HI 96701C and HI 96762C include photometer, sample cuvettes (2) with caps, 9V battery, scissors, cuvette cleaning cloth, instrument quality certificate, instruction manual and rigid carrying case.

CAL CHECK™ standards and testing reagents sold separately

REAGENTS AND STANDARDS

HI 96701

HI 96701-11 CAL CHECK™ Standard Cuvettes

HI 93701-01 Reagents for 100 tests

HI 93701-03 Reagents for 300 tests

HI 96762

HI 96762-11 CAL CHECK™ Standard Cuvettes

HI 95762-01 Reagents for 100 tests

HI 95762-03 Reagents for 300 tests

For a complete list of Reagents, see Reagents Section 18.